

Know Your Underpad

In this very competitive market, it is important to make sure you know what to look for in a quality underpad to be certain it will perform as expected.

OVERLAP.

The backing plastic should overlap the sides by AT LEAST one inch.

BACKING PLASTIC.

Should be uniform in colour. Lighter streaks in the colour or creases indicate weakness in the plastic that can lead to splitting.



CENTRE GLUE LINE.

A quality underpad should have a centre glue line to prevent the tissue slipping, keeping the underpad flat. This is easily observed on the plastic side.

SIDE GLUE LINE.

Should be thin, neat and straight, closer to the tissue than the edge to prevent the plastic folding back exposing the patient to the glue.

AUSTRALIAN MADE: Med-Con's underpads have been proudly made in Australia for 18 years. Consistent quality, service and reliability have ensured we are one of the market leaders.

PACKAGING. Underpads are best stored flat rather than folded, to prevent weakening of the backing plastic.

LATEX CONTENT. Common Latex allergies are becoming a major concern in the industry due to an increasing number of reported cases and litigations as a result thereof. Using an underpad with a latex free glue ensures the elimination of any risk associated with this product.

Phone: 03 9753 2710

Fax: 03 9753 2722

Email: info@medcon.com.au