



Med-Con Pty Ltd

Solutions for Bio-Security

Phone: 03 9753 2710

Fax: 03 9753 2722

Email: info@medcon.com.au

Sterilization Bags Storage Conditions



CSD
Dentists
Doctors
Medical
Veterinary
Tattoo Artists

To obtain the best results the products need to be stored in a cool (under 25oC), dry place out of direct sunlight.

It is recommended that stock can be used up to a 30 month period from the date of manufacture.

Stock outside of this time period needs to be monitored and tested by the end user for quality of seal. Any stock failing to seal correctly should not be used.