



## Technical data sheet

### Dust Cover Bags

Article-no.: see catalogue	Revision No.: 1
Date: 10.10.05	Page: 1 of 1

<b>Field of application:</b>	Dust Cover bags are destined for transport and storage of medical devices in their primary or secondary packaging. Note: Sterile packs must be cool and dry before packed into dust cover bags.
<b>Material:</b>	Polyethylene, transparent, 60 µm
<b>Print:</b>	Plain blue, Pantone 286
<b>Adhesive:</b>	Synthetic rubber, doublesided, on flap. The self seal closure is ten times reusable.

**Technical data is submitted to the usual variations of operating procedures.  
Therefore it is not compulsory and does not exempt from validation and operating tests.**