

PRODUCT SPECIFICATION SHEET

Product Name

Surgical Face Mask with Ear Loops
Level 1

Part Number

171604



Supplier Details

Company Name: **Med-Con Pty. Ltd.**
Contact: Sales
Phone: 1800 001 435
Email: sales@medcon.com.au



Product Specifications

<i>Product Description</i>	Surgical Face Mask with Ear Loops
<i>Mask Type (eg n95, Level 1, 2 or 3)</i>	Level 1
<i>HPV Product Number (if contracted Product)</i>	N/A
<i>Country of Origin</i>	Australia
<i>Applicable Directive & Product Standards (UNSPSC Classification)</i>	UNSPC No. 42131713 AS4381:2015
<i>(ARTG Number)</i>	ARTG No. 109872
<i>Material</i>	Polypropylene / Polyethylene Nonwoven
<i>Ply</i>	3
<i>Type of Attachment (e.g tie or ear loop)</i>	Elastic Ear Loop
<i>Latex Content</i>	Nil
<i>Visor (e.g with or without)</i>	Without
<i>Bacterial Filtration Efficiency (BFE)</i>	=>98%
<i>Fluid Resistance</i>	80mmHg
<i>Differential Pressure (Delta P)</i>	>4.0 h ₂ O
<i>Test Reports</i>	Nelson Laboratories

PRODUCT SPECIFICATION SHEET



<i>Quality Assurance</i>	ISO 9001 – Certificate No. 191010QMS – Issued Oct 19
<i>Measurements</i>	150mm x 180mm
<i>Colour</i>	White Inner – Green Outer
<i>Single Use (Y/N)</i>	Y – Single use only do not re-use.
<i>Donning and Doffing</i>	As per hospital guidelines
<i>Use</i>	If mask becomes damp or is damaged dispose of immediately If irritation occurs cease use immediately As per employers PPE Policy / Procedure.
<i>Packing (e.g Pieces/Pack, Pack/Carton)</i>	50 pieces per inner – 6 inners per outer carton.
<i>Shelf life for e.g 12 months from manufacture date</i>	5 Years from Date of Manufacture if stored in a cool, dry & clean environment. Storage temperature between -5°C - 50°C. Below 65% humidity.
<i>Disposal stream (general / infectious)</i>	General